

Resinex AP UB OH

Strong base anion exchange resin

ResinexTM AP UB OH is a high purity, premium grade, strongly basic macroporous anion exchange resin type 1. The macroporous crosslinked matrix offers a very high resistance to physical breakage and organic fouling. Its remarkable physical stability makes it highly suitable for industrial applications at very high velocities, such as condensate treatment and reversible removal of organics.

The selected bead distribution of Resinex™ AP UB OH - very close to monodisperse - is especially adapted for all modern counter-current systems (i.e. Schwebebett, UPCORE,..) and mixed bed systems.

Typical Properties

Туре	Crosslinked polystyrene divinylbenzene
Form	macroporous, milky white, spherical beads
Functional group	Quaternary amine, Type 1
Whole bead count	95% min.
lonic form, as shipped	OH-
Bead size	(≥ 90%) 0.50 - 0.71 mm
Uniformity coefficient	1.20 max.
Bulk density, as shipped	660 kg/m³
Real density	1.08 g/cm ³
Water retention	50 - 60%
Total capacity (Cl- form)	1.15 eq/l min.
Volume change OH ⁻ -> Cl ⁻	-20% max.
Stability, temperature	60°C (OH: form) max.
Stability, pH	0 - 14

Standard Design Conditions

Bed depth	> 700 mm
Service flow rate	8 - 55 BV/h
Backwash expansion	50 - 75%

Key Features and Benefits

- High Integrity Beads
 Excellent resistance to mechanical degradation ensures low pressure drop
- Excellent Resistance To Organic Fouling Removable organics
- Resistance To Osmotic Shock
 Extended lifetime and very low number of broken beads
- Uniform Bead Size
 Lower pressure drop and regenerant consumption

Typical Applications

- Demineralisation in industrial water treatment systems, especially in the presence of high organic loadings
- Demineralisation and polishing when used in combination with Resinex™ KP
- Treatment of electroplating rinse waters in combination with Resinex™ KP

Standard Packaging


- 25 lit. PE valve bag
- 1000 litre big bag






Resinex™ AP UB OH Strong base anion exchange resin

Pressure Drop



Backwash Expansion



Counter-Flow

Standard Regeneration Parameters

Concentration	4% NaOH	2% NaOH
Level	60-150 g/l	50-80 g/l
Flow rate regenerant	4-6 BV/h	6-8 BV/h
Contact time regenerant	30-60 min.	20-40 min.
Flow rate slow rinse	4-6 BV/h	6-8 BV/h
Slow rinse water required	2-4 BV	2 BV
Flow rate fast rinse	10-30 BV/h	10-30 BV/h
Fast rinse water required	6-10 BV	6-10 BV

Co-Flow

The use of a weak base solution such as ammonia or sodium carbonate as a regenerant is an alternative to caustic soda. Please contact your nearest Jacobi Carbons sales office for further information.

Product Packing



25 lit. polyethylene valve bag 48 bags per pallet



Polypropylene FIBCs (big bag), 1.000 lit.



NOTICE Jacobi Carbons reserves the right to change product specifications without prior notification. The information contained in this datasheet is intended to assist a customer in the evaluation and selection of products supplied by Jacobi Carbons. The customer is responsible for determining whether products and the information contained in this document are appropriate for the customers use. Jacobi Carbons assumes no obligation or liability for the usage of the information in this datasheet, no guarantees or warranties expressed or implied, are provided. Jacobi Carbons disclaims responsibility and the user must accept full responsibility for members based on this data.

CAUTION Strong oxidizing agents such as nitric acid can react violently with ion exchange resins and cause explosive type reactions. Before using strong oxidants, consult sources knowledgeable in the handling of these materials.





info@lenntech.com Tel. +31-152-610-900 www.lenntech.com Fax. +31-152-616-289

